

Stim Care

Innovative portable electrostimulator for aesthetics

Esthetic current portable device Stim Care

Stim Care is the most complete portable device in the market, that allows many esthetical treatments by biphasic and polarized currents: **TENS (Transcutaneous electrical nerve stimulation)**, **FES (Functional Electrical Stimulation)**, **High Force**, **Russian**, **Lipolysis**, **Galvanic**, **Micro galvanic** and **MENS**.

Main treatments:

- Pain control
- Cellulitis
- Fat reduction
- Cutaneous permeation
- Lymphatic drainage
- Rehabilitation
- Muscle sagging
- Tissue repair
- Wrinkles
- Stretch marks



Stim Hair

Platform for hair treatment



Trichology device Stim Hair

Stim Hair is a trichology device, specific for hair treatment, with the resources: **Phototherapy (LED), High frequency and Vaccum Therapy**

Main treatments:

- Hair Therapy
- Photodynamic Therapy
- Scalp cleaning
- Seborrhoeic Dermatitis

Fluence MAXX


High technology in LED and LASER multitherapy




FLUENCE MAXX is a phototherapy device with the highest LED and LASER technology. It presents a wide range of resources for multitherapy and exclusive accessories, with 28 options: Laser Pens with greater power, LED cluster (circular and linear), LED cluster + Laser, MAXX LED cluster, LED Cluster + Laser MAXX, Facial LED Phototherapy Mask, Facial LED Phototherapy Mask + Neck and Bracelet for ILIB Technique. It has more of 1,700 * pre-programmed protocols and the possibility of recording new and an INFO menu with a tutorial for to use and care for better use of all the equipment's resources.

Main treatments:

- Acne
- Acupuncture
- Alopecia
- Analgesia
- Cellulitis
- Fat Located
- Hydration
- Lymphatic Drainage
- Lymphedema
- Melasma
- Photodynamic Therapy
- Physical rehabilitation
- Rejuvenation
- Tissue repair
- Transcutaneous ILIB
- Rejuvenation



Sonic Compact MAXX *Ultrasound and Current – Body and Facial*



Sonic Compact MAXX is an innovative intelligent system for combination therapy with Ultrasound and Currents for facial and body use. This device presents the first 5MHz facial Ultrasound applicator in Brazil. The most complete and versatile of the market. At the frequency of 3 MHz, it presents intelligent calculation programming of the application time and automatic configuration of the ultrasonic intensity according to the fat layer. In addition, Sonic Compact MAXX has 2 output channels for therapeutic currents, so the equipment allows complete treatments by the resources: **1 and 3 MHz body ultrasound, Facial ultrasound, TENS, FES, HighForce, Russian, MediumForce, SoftForce, Lipolysis, Interferential 2k and 4k, Iontophoresis (Galvanic current), Microgalvanic and MENS.**

Main treatments:


- Pain control
- Cellulitis
- Fat reduction
- Cutaneous permeation
- Muscle toning
- Rehabilitation
- Stretch marks
- Tissue repair
- Water retention (edema)
- Rejuvenation



SCAN ME

Stimulus Face MAXX

Facial multi-therapies device



Stimulus Face Maxx is a facial multi-therapies device, the most complete and versatile multiplatform in the market. With all the features of the exclusive HTM Match* app, it is possible to synchronize the Stimulus Face Maxx via Bluetooth to your cellphone, which will ensure greater customization in service, assistance and updates regarding the equipment, in addition to all the other advantages that the HTM Match makes it possible. *Feature available only for Brazilian territory

The device has 4 output channels and different technologies that allow complete treatments by the resources: **HighForce current, Russian current, Galvanic current, Microgalvanic current, MENS current, Vaccumtherapy, Ultrasonic peeling, High frequency and Phototherapy**

Main treatments:

- Acne
- Cutaneous permeation
- Facial drainage
- Muscle toning
- Photodynamic Therapy
- Skin cleansing
- Tissue repair
- Rejuvenation



Empower

Technology in favor of the body and beauty



EMPOWER is a device that has High Intensity Electromagnetic Field technology. This technology is intended for muscle strengthening and body definition, capable of performing muscle contractions at maximum intensity. These maximum contractions create a high demand for energy supply to the muscles, activating the lipolytic cascade, consequently reducing the fat deposits in the applied area. EMPOWER has 4 options of ready treatments divided into beginner and advanced that allow you to achieve different types of indications: **HIIT, Performance, Phases Mode for strength or definition and Cycles Mode for strength or endurance.**

Main Treatments:

- Fat reduction
- Muscle toning
- Rehabilitation
- Pelvic Floor Dysfunction
- Urinary Incontinence



Body focused ultrasound – Body piezoelectric shock waves and Stationary and dynamic microfocused facial ultrasound



Platform for Face Lifting and Body Lipomodeling

The new ULTRAFOCUS® is a High-Intensity Focused Ultrasound (HIFU) device, developed for non-invasive aesthetic treatments. It has 2 applicators: 1 facial applicator with 3 tips and 1 body applicator with 4 exclusive tips that allow depth variation of the focal point, thus serving a wide range of applications, with automatic parameter adjustment according to the selection of the tip and application technique. The new ULTRAFOCUS® has an 8" Touch Screen display with interactive navigation and smart treatment area selection. Includes body and facial demarcation guides, making parameters use and adjustment easier.

Main treatments:

- Localized Fat
- Cellulitis
- Facial Aging
- Stretch marks
- Skin laxity (skin tightening)



Ácrus

Advanced Triple Wave LASER technology for epilation of all skin phototypes



ÁCRUS uses laser technology for hair removal with a Triple Wave applicator. The applicator allows simultaneous emission of 3 wavelengths lasers: Alexandrite (755 nm), Diode (808 nm), and ND-Yag (1064 nm) being suitable for all skin phototypes, including tanned ones. Also adapting to wavelengths for all hair thicknesses and densities. With a 12-mm² spot, providing greater treatment effectiveness and speed, the applicator has smart cooling – Cryocooling MAXX up to -10 degrees C – which ensures more comfort and safety to the patient. ÁCRUS has a 10.4' touch screen display with all the interactivity of the Pro-Performance system. ÁCRUS allows the operator to provide treatment in three different ways, namely: scanning, sequential, and spot mode, thus ensuring protocol flexibility and application speed. Equipped with all the features of the HTM Match application, ÁCRUS is an interactive device, with close and customized communication with the factory. With a unique i-Check equipment integrity system, it is ready to be used safely and efficiently.

Main Treatments:

- Hair removal





HTM
TECNOLOGIA DE RESULTADO

TRAINING VIDEOS WITH ENGLISH SUBTITLES

Ácrus: https://youtu.be/KVa0L6np_k8

Ultrafocus: <https://youtu.be/9ZCF0aKpAdw>

Sonic Compact MAXX: <https://youtu.be/-atDX7gUGyo>

Stimulus Face MAXX: <https://youtu.be/oSEubmcXfWs>

Fluence MAXX: <https://youtu.be/k2qnaNWj9Eo>



Empower  **RESULTS**

3 SESSIONS – 3 times a week

BEFORE



AFTER



BEFORE



AFTER



RESULT MADE BY PRISCILA ARAUJO

3 SESSIONS – 3 times a week

BEFORE



AFTER

BEFORE



AFTER



RESULT MADE BY PRISCILA ARAUJO

8 SESSIONS – twice a week


BEFORE

AFTER



RESULT MADE BY FLAVIA MEDEIROS

8 SESSIONS – twice a week





UltraFOCUS RESULTS

BEFORE



AFTER 3 months



4 SESSIONS - once a week
Tips: 10 mm – 15 mm – 20 mm

BEFORE





AFTER



4 SESSIONS – once a week
Tips: 10 mm – 15 mm – 20 mm

5 sessions

8 sessions





BEFORE



AFTER



BEFORE



AFTER



A: before the procedure

B: Immediately after the procedure

Tips: 4,5 mm – 3,0 mm – 1,5 mm

BEFORE



AFTER



A: before the procedure

B: Immediately after the procedure

Tips: 3,0 mm – 1,5 mm

3 sessions

BEFORE

AFTER



BEFORE

AFTER



Dynamic microfocused facial ultrasound


1 session

Tips: 4,5 mm – 3,0 mm – 1,5 mm



Ácrus  **RESULTS**

2 SESSIONS – once a month



2 SESSIONS – once a month



5 SESSIONS – once a month



3 SESSIONS – once a month



BEFORE



AFTER

Hibridi RESULTS

Same technology as Sonic Compact Maxx*

4 SESSIONS – once a week


Protocol: Lipo 3D + Stimulus



BEFORE



AFTER 3 SESSIONS



AFTER 4 SESSIONS

4 SESSIONS – once a week

Protocol: Lipo 3D + Stimulus



BEFORE



AFTER 4 SESSIONS



Fluence *MAXX* **RESULTS**

BLUE LED PROTOCOL



BLUE LED AND RED LASER PROTOCOL





Stim Hair **RESULTS**

RED LED AND VACCUM THERAPY PROTOCOL



BEFORE



AFTER