











Device	Basic spec	Treatments areas	Fields of application	Link	Picture	Retail price
Stimulus Face MAXX with all accessories + 2 Led cluster (red and blue)	<ul style="list-style-type: none"> • Biphasic currents: HighForce and Russian (120mA) - 4 channels; • Polarized currents: Galvanic (30mA), Microgalvanic (1000 µA) e MENS (1000µA) - 2 channels; • Vacuum (-70kPa/525mmHg); • Ultrasonic peeling (26kHz); • High frequency; • Phototherapy; • 87 pre-programmed protocols; • Higher number of phototherapy applicators (red, blue, infrared, green, amber, violet). 	Face, Neck, Hair, Body (some specific regions).	Aesthetic; Trichology (hair); Podiatry; Pre and post micropigmentation/ablative laser / mechanical peels.	https://htmeletronica.com.br/en/products/stimulus-face-maxx-facial-multi-therapies-device/		\$ 3,450.00
Stim Hair	<ul style="list-style-type: none"> • Phototherapy (red + infrared); • Vacuum (-70kPa/525mmHg) with specific accessories for hair treatment; • High frequency. 	Face, Neck, Hair.	Trichology (hair).	https://htmeletronica.com.br/en/products/stim-hair-trichology-device/		\$ 1,800.00
Beauty Steam MAXX	<ul style="list-style-type: none"> • Operates on two mains voltages - bivolt (110V~/220V~) - first in the world; • Features cluster of LEDs for Chromotherapy and Photobiomodulation, allowing the generation of different colors (Blue, Green, Red, Violet, Yellow and Amber); • More than 50 pre-programmed protocols and possibility to record new ones; • Compartment for essential oils. 	Face, Neck, Hair.	Aesthetic; Chromotherapy Trichology (hair).	https://htmeletronica.com.br/en/products/beauty-steam-maxx-3/		\$ 1,140.00
Sonic Compact MAXX with body and facial applicator	<ul style="list-style-type: none"> • 1 and 3 MHz body ultrasound; • Facial ultrasound - Brazil's first 5MHz facial applicator; • Biphasic currents: TENS, FES, HighForce, Russian, MediumForce, SoftForce, Lipolysis, Interferential 2k and 4k (2 channels); • Polarized currents: Iontophoresis, Microgalvanic and MENS (2 channels). 	Face, Neck, Body	Aesthetic; Physiotherapy; Rehabilitation; Sports.	https://htmeletronica.com.br/en/products/sonic-compact-maxx-3/		\$ 1,800.00

<p>Hibridi with body and facial applicator</p>	<ul style="list-style-type: none"> • 1MHz and 3MHz body ultrasound with exclusive FloatingHead applicator: perfect coupling with the applied surface; • Lightweight and ergonomic body applicator with higher output power: 63W and 21cm² ERA; • Facial ultrasound – First 5MHz facial applicator in Brazil; • Biphasic currents with greater intensity: HighForce (1000Hz), Russian (2500Hz), MediumForce (4000Hz) and SoftForce (8000Hz) – 4 channels (150mA); – Drainage – Fast, Medium and Slow Speed – 4 channels (150mA); – Transcutaneous and Percutaneous Lipolysis – 4 channels (150mA); • Polarized currents: <ul style="list-style-type: none"> – Iontophoresis – 4 channels (30mA); – MENS – 4 channels (1000uA); • Possibility of performing single therapies, combined therapies and simultaneous therapies. 	<p>Face, Neck, Body</p>	<p>Aesthetic; Physiotherapy; Rehabilitation; Sports.</p>	<p>https://htmeletronica.com.br/en/produktos/hibridi-combined-therapy-device/</p>		<p>\$ 6,300.00</p>
<p>Effect</p>	<ul style="list-style-type: none"> • Compact radiofrequency for body, facial and intimate treatments; • Rapid reach of therapeutic temperature; • Guarantee maintenance of therapeutic temperature even in large areas; • Multifrequency technology: 0.6MHz, 1.2 MHz and 2.4 MHz; • Monopolar, bipolar and multipolar application; • Resistive and capacitive tips; • Innovative Pro Performance System: designed for the user's experience. 	<p>Face, Neck, Body and intimate areas</p>	<p>Aesthetic</p>	<p>https://htmeletronica.com.br/en/produktos/effect-radiofrequency-device/</p>		<p>\$ 3,600.00</p>
<p>Beauty Dermo MAXX</p>	<ul style="list-style-type: none"> • Allows you to choose between treatments: Vacuum LED, Vacuum Therapy, Phototherapy, Pressotherapy, Cupping Therapy, Diamond Peeling and Buttock Pumping; • Features Vacuum LED applicator containing: 3 green LEDs (300mW), 3 red LEDs (300mW) and 3 infrared LEDs (300mW), in continuous and pulsed modes. Possibility of specific or combined use of LEDs, with or without the use of the vacuum; • New key for Cupping Therapy, with special and selective suction system for application of suction cups; • Vacuum operating mode: Continuous and Pulsed (5 to 50 pulses/min); • New digital pressure sensor with variable adjustment of negative pressure from -1.5kPa to -80kPa* (*equivalent to -600mmHg); • Positive pressure adjustment for Pressotherapy up to 250mmHg; • Pressotherapy with 2 output channels: CHANNEL A for legs or arms and CHANNEL B for abdomen. 	<p>Face, Body</p>	<p>Aesthetic; Rehabilitation; Sports.</p>	<p>https://htmeletronica.com.br/en/produktos/beauty-dermo-maxx-3/</p>		<p>\$ 4,200.00</p>

<p>Fluence MAXX (with mask, 1 laser probe, 1 cluster maxx (laser + LED))</p>	<ul style="list-style-type: none"> • High technology in LED and LASER multitherapy; • Individual or associated therapies; • Intelligent programming with directed protocols; • Different models of masks for phototherapy; • New MAXX Clusters options with larger application area; • Exclusive a bracelet for ILIB technique. 	<p>Face, Body</p>	<p>Aesthetic; Trichology (hair); Podiatry; Pre and post micropigmentation / ablative laser / mechanical peels; Physiotherapy; Rehabilitation; Sports.</p>	<p>https://htmeletronica.com.br/en/producos/fluence-maxx-3/</p>		<p>\$ 4,650.00</p>
<p>Stimulus Physio MAXX</p>	<ul style="list-style-type: none"> • Biphasic Currents: TENS, FES, HighForce, Russian, MediumForce, SoftForce, Interferencial 2000Hz and 4000Hz; • Polarized current: Microgalvanic, MENS; • 4 independent output channels; • 51 pre-programmed protocols. 	<p>Body</p>	<p>Physiotherapy; Rehabilitation; Sports.</p>	<p>https://htmeletronica.com.br/en/producos/stimulus-physio-maxx-excitomotor-current-device/</p>		<p>\$ 1,200.00</p>
<p>Stimulus-R</p>	<ul style="list-style-type: none"> • Biphasic currents: TENS, FES, HighForce, Russian, MediumForce, SoftForce (250 mA) - 10 channels; • Polarized currents: Galvanic (30mA), Micro galvanic (1000 µA) e MENS (1000µA) - 2 channels; • 46 pre-programmed protocols. • Electroacupuncture protocol. 	<p>Body</p>	<p>Aesthetic; Physiotherapy; Rehabilitation; Sports.</p>	<p>https://htmeletronica.com.br/en/producos/stimulus-r-multicurrent-device/</p>		<p>\$ 2,100.00</p>